

# Introduction of the Imago LEAP4000X: Second-Generation Commercial Laser-Pulsed Atom Probe Tomography



Imago Scientific Instruments Corporation  
5500 Nobel Drive, Madison, WI 53711 USA

## Focused Laser Delivers Performance

The LEAP4000X builds on the success of the LEAP3000X, the first commercially available laser-pulsed atom probe (20 units sold since release in mid 2006) to offer major advances in laser mode performance.

### Key Features:

- ✓ Ultra-small, diffraction-limited spot for optimized thermal pulsing
- ✓ Significantly increased mass resolution
- ✓ New Active Laser Control™ with fully integrated optics chain to drive evaporation process and enhance laser power stability
- ✓ Advanced auto alignment protocols
- ✓ Increased dynamic range on laser power
- ✓ Increased repetition rate (up to 1 MHz)



### Benefits:

- ✓ Change-on-the-fly operation for enhanced ease of use and improved throughput
- ✓ Improved compositional accuracy and precision
- ✓ Enhanced detection sensitivity limits
- ✓ Elimination of thermal tails in low-thermal-diffusivity materials
- ✓ Enhanced specimen yield for dielectric/ceramic materials
- ⇒ Improved quantification


## Smaller Spot = Higher Quality Data

Published experimental results show that laser-induced field evaporation is overwhelmingly dominated by a thermal mechanism<sup>1,2</sup>.


An ultra-small spot design maximizes performance.



<sup>1</sup>J.H. Bunton *et al.*, "Advances in Pulsed-Laser Atom Probe: Instrument and Specimen Design for Optimum Performance," *Microscopy and Microanalysis*, Vol. 13, No. 6 (2007) pp. 418-427.  
<sup>2</sup>F Vurpillot *et al.*, "Thermal response of a field emitter subjected to ultra-fast laser illumination" *J. Phys. D: Appl. Phys.* 42 (2009) 125502 (7pp).



## Ultra-small Spot Improves Quantification



- LEAP 4000X data shows improved mass resolution.
- In this example, the Phosphorous peak (located in silicon tail) is detected with much greater sensitivity.

## Application Example: High-κ Dielectric

- Shows the benefit of superior mass resolution of the LEAP4000X for separating closely spaced peaks in a complex material
- It is also an example of quality data from a low-thermal-diffusivity material

